




Safe, simple, unique

# WinFusion



# WinFusion – Safety winged IV catheter:

**HIGHEST SAFETY WITH  
COMPLETELY CLOSED SYSTEM**



## WinFusion - Your advantage

IV access for infusion therapy, intensive care and anesthesia

### Unique system

Withdrawing the cannula into the completely closed system is easy and provides highest safety for the vascular access.

### DEHP-free



All components are free from DEHP. The position in the vessel can be made visible with the radiopaque markings.

### Integrated tubing

The direct connection between catheter and tubing with clamp increases the safety for your patient.

### Facilitated Cannulation

The integrated wing with maximal grip facilitates cannulation and fixation. The wing can be removed by hand if required.

Diameter $\varnothing$	Length	Colour	Order Code
20 G / 1,1 mm	3/4" / 19 mm		WF20G19
22 G / 0,9 mm	3/4" / 19 mm		WF22G19
24 G / 0,7 mm	3/4" / 19 mm		WF24G19

Tubing length 35 cm, DEHP - free

**STERILE EO** WinFusion Safety IV catheters are EO-sterilized.



Bionic Medizintechnik GmbH – a JMS corporation  
 Max-Planck-Straße 21 · D-61381 Friedrichsdorf, Germany  
 Phone: +49(0)6172-75760 · Fax: +49(0)6172-757610  
 E-mail: info@bionic-jms.com · www.bionic-jms.com

