

# Bionic Accessories



AV-EXT2

Extension lines for dual lumen catheters, each 10 cm long. The clamps in red and blue assure safe and fast clamping.



SY100LL

Disposable syringe 100 ml with Luer Lock connector.



SNA3

Y-Adapter two connection tubes, each 10 cm long. The colour coded clamps assure fast and safe clamping.



CLAMP-SP-02

Stop&Protect - simple disposable compression strap, practical dispenser box with 50 pieces.



CATH-TAS

DIACATH catheter dressing with pocket made of non-woven fleece for easy and quick fixation of catheters.



K2565

Plastic scissor clamp for easy clamping of tubes of up to 7 mm diameter.



KPA06006-4 / KPA06004-4

Connection and disconnection set for single or double lumen catheter, contains different components, e.g. drape.



KPA06021-6

Connection and disconnection set for shunt puncture; contains various components, e.g. drape, plaster.



KAB75

Sharps container with a volume of 7.5 liters.




HEDclip

HEDclip for strain relief of all common blood tube systems with a tube diameter of 6 mm.



IKLL-ISO

Latex free injection cap with self-sealing injection membrane for Luer Lock connectors. Available in yellow.



UVL

Universal cap for Luer Lock connectors. Available in red.

# Bionic Accessories

Order code	Description	Packaging unit	Connector-material
AV-EXT2	Extension lines for double lumen catheters	50 units / carton	PC, PP
CATH-TAS	DIACATH Catheter dressing	30 units / carton	
CLAMP-SP-02	Stop& Protect compression strap, disposable	50 units / carton	
HEDclip	Strain-relief system HEDclip	10 units / carton	
IKLL-ISO	Injection cap LL, sterile, latex-free membrane	100 units / carton	PP
KAB75	Sharps container 7,5 l	40 units / carton	
KPA06006-4	Connection and disconnection set for single lumen-catheter	20 units / carton	
KPA06004-4	Connection and disconnection set for double lumen-catheter	15 units / carton	
KPA06021-6	Connection and disconnection set for shunts	120 units / carton	
K2565	Plastic scissor clamp	50 units / carton	
SNA3	Single-Needle-Adapter with Luer-Lock and clamp	50 units / carton	PP, PVC
SY100LL	Disposable Syringe 100 ml	25 units / carton	PP
UVL	Universal cap, Luer-Lock, red	100 units / carton	ABS

